irudek

IRUDEK 2000 S.L.
Pol. Apatta, naves C14-18
20400 Ibarra-Tolosa (Guipúzcoa) SPAIN

T (+34) 943 692 617

F (+34) 943 692 526

irudek@irudek.com

www.irudek.com

irudek
GROUP

irudek

MANUAL DE INSTRUCCIONES CAMILLA DE RESCATE
INSTRUCTION MANUAL RESCUE STRETCHER

IRUDEK ECO

CE 0120

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1. GENERAL INFORMATION

a. 1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, apart from attaining an adequate use of the appliance, will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

1.2 Conservation of the instruction manual

The instruction and maintenance manual must be kept with the product, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

2. GENERAL WARNINGS

Before carrying out any kind of operation on the product, the operator must carefully read the enclosed instructions paying attention to the correct safety precautions, the procedures to be followed for installation and correct usage.

- Regularly check the appliance. Carry out the prescribed maintenance in order to keep the appliance in good condition and to guarantee correct functioning and a long life.
- In the case of any abnormalities or damage to the appliance, which could jeopardize the functioning and safety, the appliance must be immediately removed from service.
- Do not alter or modify the appliance in any way, such interference could cause malfunctions and injury to the patient and/or rescuer.
- Ensure that all the necessary precautions are taken in order to avoid hazards that can arise as the result of contact, with blood or body fluids.
- Handle with care

2.1 Specific Warnings

- The product must be used by trained personnel only.
- The device is NOT designed as a spinal immobilization device. If a spinal injury is suspected, secure the patient to an approved spinal immobilization device prior to placing him/her in the stretcher.

NOTA: The device should be stored in the back pack, as prolonged exposure to sunlight (UV rays) can damage all plastics.

3. DESCRIPTION OF PRODUCT

3.1 Intended use

The recovery stretcher IRUDEK ECO has been configured as an extremely flexible rescue system, useable in any context where patient transport is needed and particularly useful

in morphologic and difficult climatic situations where a fast recovery is requested. This stretcher can be anchored with the help of the clips fixed to the belts, to the helicopter with the dedicated rope. The stretcher can be used as structure for lifting other equipment. The system is compatible with spinal boards, scoop stretcher, spinal supports, and many other immobilization devices.

3.2. Main components

1. Loading surface
2. Metal rings
3. Polyester rope
4. Handles
5. Lifting belts

3.3 Environmental conditions

Temperature of use: - 10° C a + 40°C

Relative humidity: -5% a +85%

3.4 Technical data

RECOVERY STRETCHER IRUDEK ECO	
Maximum load	200kg
Weight of device	5.5kg
Total weight, with bag and belts.	8kg
COMPONENTS (Loading surface)	
Length	2440mm
Width	918mm
Material	Poly Ethlene
Hooks (material)	Alloy Steel
Rope	9 meter
Length	9000mm
Diameter	10mm
Lifting belts	
Length	Feet zone 2030mm/Head zone1930mm
Width	Feet zone 5mm/Head zone5mm
Maximum load belts	3000kg
Maximum load Stitching	2000kg
Handles	
Material	Nylon
Handles on every longitudinal side	2

4. OPERATING INSTRUCTIONS FOR RECOVERY STRETCHER ECO SEKURAL

4.1 Transportation and Storage

Before transporting the appliance, make sure that it is correctly packaged ensuring that there are no risks of shocks or falls during the transport itself. Make sure the application of the supplied safety belt around the device is secure, before pulling it inside the

backpack and transporting it. Keep the original packaging for any eventual further transport. Damage to the appliance caused during transport and handling is not covered under the guarantee.

Repairs or replacements of the damaged parts are the responsibility of the client. The appliance must be stored in a dry place free from humidity.

After use, replace the stretcher inside its container, in order to prevent it from potential decays caused by elevated exposure to ultra-violet rays.

4.2 Preparation

On receiving the product:

- Remove the packing and display the material so that all is visible.
- Check that all the components/pieces on the accompanying list are present.

NOTE: The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage.

4.3. Functioning

The use of the device must take place as described hereafter:

(A) To unroll the stretcher and lay flat

- a) Remove the device from pack and place on the ground
- b) Unfasten retainer strap
- c) Step on foot end of the stretcher, and unroll completely to opposite end.
- d) Before starting the rescue operation, bend the stretcher in half and back roll. Repeat with opposite end of stretcher. The device will now lay flat.

(B) How to place a patient


1.- Method 1

- a) Place the device next to patient. Ensure head end of stretcher is adjacent to head of patient. Place cross straps under the stretcher.
- b) Place the patient and slide the stretcher under the patient as far as possible. Gently roll patient down on the stretcher.
- c) Slide patient to centre of the stretcher. Be sure to keep spinal column as straight as possible.
- d) Pull straps out from under the stretcher and fasten to buckles.

2.- Method 2

- Position foot end of the stretcher at the head of patient.
- Have one rescuer straddle the stretcher and support the patients head, neck and shoulders.
- Grasp foot straps of the stretcher and slide it under the patient.
- Lay the patient at the centre of the stretcher and fasten straps to buckles.

(C) How to position and fasten straps to buckles

- Unroll the stretcher and lay flat.
- Lift sides of the stretcher and fasten the four cross straps to buckle directly opposite the straps.
- At the end of the stretcher there will be a few unused metal rings, which can be used for fixing the bells for immobilization..
- Buckle the straps of the stretcher in proper manner.

(D) Lifting operations and positioning on the ground

Two nylon belts are used for anchoring and transport of the stretcher. The belt to be used for the upper part (head) is about 10 cm shorter than the one used for the lower Part (feet)

- Insert one end of head strap through lift slot at the head end of stretcher.
- Bring the straps under stretcher and through lift, slot on opposite side of stretcher.
- Make sure the two ends of the belt are equally long repeat the operation at the other end of the stretcher.
- Make sure the four ends of the two belts are fix them to connector.

E) Vertical lift

Use the supplied rope length of 9 meter with figure of 8 knot in the centre.

- Pass each end of the rope through grommets at the head of the stretcher. Pull the knot up against the stretcher.
- Continue feeding rope through unused grommets and carrying handles all the way to the foot end of stretcher. Ensure both ends of the rope are even.
- Pass the rope ends through grommets at the foot end of stretcher from the inside outward. Tie the ends of the rope together with a square knot.
- Bring ends of the rope up over the end of the stretcher. Pass through carrying handled and secure with a square knot safely each side with an over hand knot.

IMPORTANT: *The device, useable for rescue operations in narrow spaces, heights, rough terrains and in water, may never be lifted with one ring only.*

MANUAL TRANSPORT: With the help of the set of four handles it is possible to transport the device manually. The free metal rings of the stretcher can be used as anchorage for the handles. Transport belts can be connected to the handles in order to side the stretcher on the terrain.

ROLLING UP THE DEVICE

Lay the stretcher out and place retainer strap, buckle side down under foot end of the stretcher. Starting at the head end, roll the stretcher up as tight as possible. Continue to roll the stretcher up using knee to keep stretcher from unrolling. Fasten pre-placed retainer strap to buckle and place the stretcher in the back pack.

PROBLEMS	CAUSE	REMEDY
Damages to the loading surface	Improper use	Take the device out of service immediately and contact the service centre.
Damages to the metal rings dedicated the various belts and ropes	Improper use	After every use, check the conditions of the to rings. In case permanent deformations are detected, take the device out of service immediately and contact the sercie centre.
Damages to belts and ropes	Prolonged use	Periodically change components dedicated to extrication and lifting.

5.MANTENIMIENTO

5.1 Cleaning

The loading surface as well as belts and ropes, may be contaminated during various phases of the rescue and transport operations, by body liquids and blood. With the help of the materials used for the constructions of the device like polymers and other synthetic materials, it is possible to use, without any risk of incompatibility, usual products for disinfection. Dis-infection should be carried out with normal disinfectants for example any anti-septic lotion, every time after use.

5.2 Maintenance

- The device does not require scheduled servicing.
- Before and after use, in order to detect any anomalies or damages, verify the correct functioning of the device by identifying and removing anomalies which could compromise its integrity and correct operation.

CHECK CARD

The check card should be filled in before the first use.

All information about the personal protective equipment (name, serial number, date of purchase and date of putting into operation, user name, periodic examination and repair history, and next periodic examination date) shall be noted into de check card.

All records in the check card can only be filled in by a competent person.

Do not use the personal protective equipment without a duly filled check card.

EQUIPMENT RECORD

REFERENCE	<input type="text"/>
LOT, SERIAL NUMBER	<input type="text"/>
YEAR OF MANUFACTURE	<input type="text"/>
IT CAN BE USED WITH	<input type="text"/>
DATE OF PURCHASE	<input type="text"/>
DATE OF FIRST PUT INTO USE	<input type="text"/>
USER NAME	<input type="text"/>
COMMENTS	<input type="text"/>

PERIODIC EXAMINATIONS AND REPAIR HISTORY

DATE	REASON FOR ENTRY (Periodic examination/Repair)	COMPETENT PERSON Name + Signature	COMMENTS	NEXT DUE DATE FOR Periodic Examination